

Module 4

Considerations for Parenteral Products

ICH Q3D Elemental Impurities

Disclaimer:

This presentation includes the authors' views on Elemental Impurities theory and practice. The presentation does not represent official guidance or policy of authorities or industry.

International Council for Harmonisation of Technical Requirements for
Pharmaceuticals for Human Use

Legal Notice

This presentation is protected by copyright and may be used, reproduced, incorporated into other works, adapted, modified, translated or distributed under a public license provided that ICH's copyright in the presentation is acknowledged at all times. In case of any adaption, modification or translation of the presentation, reasonable steps must be taken to clearly label, demarcate or otherwise identify that changes were made to or based on the original presentation. Any impression that the adaption, modification or translation of the original presentation is endorsed or sponsored by the ICH must be avoided.

The presentation is provided "as is" without warranty of any kind. In no event shall the ICH or the authors of the original presentation be liable for any claim, damages or other liability arising from the use of the presentation.

The above-mentioned permissions do not apply to content supplied by third parties. Therefore, for documents where the copyright vests in a third party, permission for reproduction must be obtained from this copyright holder.

Content

- **Wording in guideline**
- **Principles supporting the 2 L parenteral limit**
- **Examples**

Q3D: Large Volume Parenteral Section

Parenteral drug products with maximum daily volumes up to 2 liters may use the maximum daily volume to calculate permissible concentrations from PDEs. For products whose daily volumes, as specified by labeling and/or established by clinical practice, may exceed 2 liters (e.g., saline, dextrose, total parenteral nutrition, solutions for irrigation), a 2-liter volume may be used to calculate permissible concentrations from PDEs. (*Reference 4*)

4) Holliday MA, Segar WE. The maintenance need for water in parenteral fluid therapy. *Pediatrics* 1957;19:823-32.


Maintenance Fluid Replacement

- In 1957, Holliday and Segar proposed a method for determining the average daily fluid replacement need for hospitalized patients. (Holliday, Segar [1957], ref. 1)
 - Summary, based on patient body weight
 - First 10 kg body weight: 100 ml/kg/day
 - Second 10 kg body weight: 50 ml/kg/day
 - Thereafter: 20 ml/kg/day
- This method is still considered the standard of care.
 - “The Holliday-Segar equation remains the standard method for calculating maintenance fluid requirements.” (Meyers, [2009], ref. 2)

1. Holliday MA, Segar WE. The maintenance need for water in parenteral fluid therapy. *Pediatrics* 1957;19:823-32.
2. Meyers, RS. Pediatric Fluid and Electrolyte Therapy, *Journal of Pediatric Pharmacology Therapeutics* 2009;14:204-11

Holliday-Segar nomogram

- Summary: 1 ml of replacement fluid is required for every 1 Calorie expended. Caloric expenditure is correlated with body weight.
- A 50 kg patient requires about 2000 ml of fluid per day parenterally.
- Q3D established 2 liters as the reference volume for determining elemental impurity levels in parenteral products that may be delivered in large volume.
 - Q3D PDEs are based on 50 kg body weight
- On the basis of body weight, the Holliday-Segar nomogram confers an additional safety factor when volumes exceeding 2 liters are delivered to patients whose body weight exceeds 50 kg.



Example 1: 10% Dextrose Injection, USP

- **Product information**
 - Supplied in 1L PVC bags
 - Dose1: daily dose set at 2L
 - USP water for injection is used to produce the drug product
- **Risk assessment approach**
 - Risk assessment conducted using the drug product approach
 - During a review of the components of the drug product no potential elemental impurities were identified as likely to be present
 - Manufacturer conducted testing on samples from 3 representative commercial batches of the drug product to support the risk assessment (testing included class 1, 2A and 3 elements)
 - Leachable studies demonstrate that there are no significant elemental impurity contributions from the container closure system over the shelf life of the product.

1. Package Insert Dosage: “As directed by a physician. Dosage is dependent upon the age, weight and clinical condition of the patient.”

Example 1: 10% Dextrose Injection, USP

Lot number	As, µg/L	Cd, µg/L	Hg, µg/L	Pb, µg/L
20xx001	<0.50	<0.10	<0.25	<0.10
20xx002	<0.50	<0.10	<0.25	<0.10
20xx003	<0.50	<0.10	<0.25	<0.10
	As, µg	Cd, µg	Hg, µg	Pb, µg
Total daily level ¹	<1.0	<0.20	<0.50	<0.20
PDE	15 µg/day	2 µg/day	3 µg/day	5 µg/day
Control threshold	4.5 µg/day	0.6 µg/day	0.9 µg/day	1.5 µg/day

Lot number	Co, µg/L	V, µg/L	Ni, µg/L	Li, µg/L	Sb, µg/L	Ba, µg/L	Mo, µg/L	Cu, µg/L	Sn, µg/L	Cr, µg/L
20xx001	<0.20	<0.40	<2.0	<10.0	<2.0	<1.0	<10.0	<2.0	<3.0	<10.0
20xx002	<0.20	<0.40	<2.0	<10.0	<2.0	<1.0	<10.0	<2.0	<3.0	<10.0
20xx003	<0.20	<0.40	<2.0	<10.0	<2.0	<1.0	<10.0	<2.0	<3.0	<10.0
	Co, µg	V, µg	Ni, µg	Li, µg	Sb, µg	Ba, µg	Mo, µg	Cu, µg	Sn, µg	Cr, µg
Total daily level ¹	<0.40	<0.80	<4.0	<20.0	<4.0	<2.0	<20.0	<2.0	<6.0	<20.0
PDE	5 µg/day	10 µg/day	20 µg/day	250 µg/day	90 µg/day	700 µg/day	1500 µg/day	300 µg/day	600 µg/day	1100 µg/day
Control threshold	1.5 µg/day	3 µg/day	6 µg/day	75 µg/day	27 µg/day	210 µg/day	450 µg/day	90 µg/day	180 µg/day	330 µg/day

¹ (2L x observed level)

Example 1: 10% Dextrose Injection, USP

- **Conclusions**

- All values for the potential elemental impurities of concern were observed to be less than the detection limit.
- For each element, the limit of quantitation is less than the control threshold.

Example 2: 0.45% Sodium Chloride Injection, USP

- **Product information**
 - Supplied in 1L PVC bags
 - Dose¹: daily dose set at 2L
 - USP water for injection is used to produce the drug product
- **Risk assessment approach – component approach**
 - An assessment of the water for injection revealed that the only potential elemental impurities to be considered were As, Cd, Hg, and Pb.
 - Literature revealed mean results for three grades of salt using validated methods (<http://eusalt.com/salt-quality>, As, Cd, Hg, Pb testing in salt)
 - Manufacturing equipment does not contribute elemental impurities due to mild conditions employed in the manufacture of the drug product
 - Leachable studies demonstrate that there are no significant elemental impurity contributions from the container closure system
 - Monitoring of the WFI quality demonstrates that there are no significant elemental impurity contributions.

¹ Package Insert Dosage: “As directed by a physician. Dosage is dependent upon the age, weight and clinical condition of the patient.”

Example 2: 0.45% Sodium Chloride Injection, USP

	As, conc.	Cd, conc.	Hg, conc.	Pb, conc.
Maximum level in NaCl*	0.024 µg/g	0.011 µg/g	0.00062 µg/g	0.0808, µg/g
Water for injection	<0.50 µg/L	<0.10 µg/L	<0.25 µg/L	<0.10 µg/L
Impurity level per liter of drug product	<0.5108	<0.1495	<0.25279	<0.464
	As, µg	Cd, µg	Hg, µg	Pb, µg
Total µg/day (based on total 2L/day)	<1.0216 µg	<0.299 µg	<0.50558 µg	<0.928 µg
PDE	15 µg/day	2 µg/day	3 µg/day	5 µg/day
Control threshold	4.5 µg/day	0.6 µg/day	0.9 µg/day	1.5 µg/day

* NaCl data reported represent the maximum levels reported in the literature (see references)

Conclusions:

- **All potential elemental impurity levels are lower than the control threshold for each elemental impurity**

References

- **European Committee for the Study of Salt, ECSS/CN 287-1982, Statistical Study of Inter-Laboratory Analysis of Sodium Chloride (As, Cd, Hg, Pb)**
 - As: EuSalt/AS 011-2005
 - Hg: EuSalt/AS 012-2005
 - Pb: EuSalt/AS 013-2005
 - Cd: EuSalt/AS 014-2005